



## **UNISEX Compression Socks with BAMBOO SmarTec® fabric**

*Made of Bamboo, Nylon, cotton and elasthan, BAMBOO SmarTec® fabric is designed for people who demand accurate gradual compression socks with extra skin care solutions*

*BS 6612:1985 CERTIFIED GRADUAL COMPRESSION 25mmHg  
100% ORGANIC ORIGIN BAMBOO FIBRE*

**Certified Gradual Compression Socks effects and benefits**



Compression socks avoid ankle swelling and are indicated as a valid varicose veins treatment for:

- Stage I (telangiectasia or reticular veins)
- Stage II (varicose veins with no symptoms or moderate symptoms)
- Some skin changes due to venous disorders, such as eczema/dermatitis, specifically around the ankle area

Gradual compression can properly improve blood flow and muscle oxigenation, pushing the returning blood from the bottom of your legs towards your heart.

**Bamboo fibre Skin Care**

Bamboo SmarTec® composition is 50% organic bamboo fiber, 10%Nylon, 20% Cotton and 20% Elastan.

Bamboo fiber has not been chemically treated so to preserve all its properties:

- extremely soft to the touch: reduces skin irritation
- highly absorbent: regulates skin temperature
- natural antibacterial: reduces perspiration, reduces allergical skin symptoms



Elastan is added to obtain the compression effect and Nylon and Cotton to reinforce structure.

**UNISEX available colors and sizes**



black



blue sky



orange

INDESmed compression socks are produced with 25mmHg (100% compression) in the ankle area, 17.5mmHg (70% compression) in the calf, and 12.5mmHg (50% compression) in the popliteal fossa under British Standard 6612:1985 directive.

|          |  |          |   |          |  |
|----------|--|----------|---|----------|--|
| <b>L</b> | US WO'S: 12,5-16<br>US MENS: 10,5-14<br>UK: 10-13,5<br>foot measure: 10,75"-11,87" | <b>M</b> | US WO'S: 8,5-11,5<br>US MENS: 7-10<br>UK: 6,5-9,5<br>foot measure: 9,62"-10,56" | <b>S</b> | US WO'S: 5-8<br>US MENS: X-6,5<br>UK: 3-6<br>foot measure: 8,50"-9,50" |
|----------|--|----------|---|----------|--|

**Compression socks dressing and washing instructions**

Put on your compression socks starting by sliding your toes into the rolled sock and then pull the botton over your heel.

cold water 80°F/ 30°C – gentle cycle

