

Crutches for Kids



Aluminum Colored Forearm Crutches for Kids

INDESmed crutches for kids. Lightweight adjustable crutches.

INDESmed crutches for kids keep the natural position of the child's wrist avoiding tendinitis and hand and arms pain.

Curved and oval shape allows a perfect regulation and anatomic adaptation.

Continuous regulation of the handgrip allows accurate adjusting according to the height of your child increasing comfort.

Handle silicone covers are made of soft medical grade silicone (FDA certified).

These crutches are silent, with no clicks when walking.

Made of sturdy aluminum, tested according to EEC norms, INDESmed crutches for kids are available in multiple amazing colors.

They are adaptable to the millimeter to the height of the kid.

Crutches' sizes for kids

How do I choose the correct sizing for my kid crutches?

Every child's forearm crutch is fully adjustable simply moving the cuff and handgrip along the curved area.

Technical data

- Crutch weight: 0.89lb
- Certification: CE 11344-1.2007 - FDA LISTED
- Max load: 176lb (80kg) / crutch



	Weight	Material	Certificates	Features
Crutch Structure	221gr	Aluminium	CE 100Kg load test	Oval shape
Crutch Handle	90gr	Nylon	CE 200Kg load test	Height adjustable
Handle Cover	10gr	Silicone	FDA approved	4 couleurs, interchangeables
Crutch Open Cuff	90gr	Nylon	CE	Height adjustable d1 (7cm KIDS) d2 (8 cm STANDARD) d3 (9 cm BIG)
Crutch Closed Cuff	120gr	Nylon	CE	Height adjustable
Crutch tip	25gr	High Friction Rubber	CE 200Kg load test	22° angle



Designed and manufactured by VIDCAPP ITHEF SL (INDESmed®) in Madrid – Spain

Technical Service: info@indesmed.com

www.indesmed.com