## INDESmed Improving the Future



## Ultralight ergonomic crutches designed for style, engineered for freedom

INDESmed patented design offers an accurate body adaptation system to match your everyday crutch-aid requirements. The unique curved shape improves comfort, avoiding tendinitis and other injuries while the ergonomic padded handgrip helps preventing calluses on your hands.

Our crutches are available in four colors. Get them! You won't stop moving!



User's height	Recommended size
175-185 cm	L (h = 110 cm)
165-175 cm	M (h = 105 cm)
155-165 cm	S (h = 100 cm)

EPDM and Nylon-66 tip

Three different cuff sizes (2.75", 3.14" or 3.54" in diameter)

## A crutch that perfectly adapts to you

Choose your crutch size in the box above and feel how its exclusive curved style gets your exact body measurements while its fixing system avoids the annoying crutching noise.

Besides, there are four fantastic aluminum colors combining with those for the grip cover!

## No more tendinitis

Keeping your wrist's natural position, the gripping angle assures you will not suffer tendinitis or any other pain on your hand, such as medium or ulnar nerves pain.

The handgrip spreads your weight in an optimal way on your hand, so you will not have calluses anymore after a long period using your crutches. The hypoallergenic silicone cover is available in four colors and it is incredibly soft.

A high friction coefficient EPDM tip, together with its high contact angle, increases the user's safety.

Choose the silicone color

www.indesmed.com



EN-ISO 11344-1:2007