

Ultralight ergonomic cuff

Technical specifications

| | | |
|---|----|------------------------------|
| Weight | | 25 grams |
| Material | | black-coloured ABS. |
| Internal diameter <i>Three sizes available</i> | d1 | 7cm/2.75" – Kids. |
| | d2 | 8cm/3.14" – Standard, adult. |
| | d3 | 9 cm/3.54" – Large. |



1. Ultra-light ergonomic cuff, made in black ABS.
2. Its reduced weight helps to the easy use of the crutch.
3. Continuous height regulation along the curvature of the rincipal structure of the crutch.



figure 1

4. The diameter of the cuff, available in different sizes, improves the user's adaptation. See figure 1.