

A toucan bird with a large, colorful beak (yellow, orange, and red) and black body is perched on a wooden branch. To its right is a walking cane with a black carbon fiber shaft and a bright orange handle. The background is a light, out-of-focus scene with more branches.

Ultralight, Carbon fiber
walking cane

A new bird has arrived
to the jungle!

Technical specifications

Total weight	150g / 5.2 oz (\pm 5%)
Maximum load	100Kg / 220lb
Sizes	There are four sizes available according user's height.
Main frame material	100% Carbon fiber.
Load cycles	1.10^6 (Minimum)
Gripping angle (β)	($\beta=104^\circ$)

A modern and fashionable walking cane with an ergonomic handgrip and a soft silicone cover.

It's exclusive curved and oval carbon fiber frame allows optimal load transfer, and it's ergonomic handgrip ensures proper hand positioning avoiding wrist pain.



Sizes and colors description

Advisable sizes according user's height		Cane, Total length
Size XL	From 73" (185cm) to 77" (195 cm)	h=98cm(46")
Taille L	From 69" (175cm) to 73" (185 cm)	h=88cm(43")
Taille M	From 65" (165cm) to 69" (175 cm)	h=83cm(41")
Taille S	From 61" (155cm) to 65" (165 cm)	h=78cm(39")

Soft Silicone cover colors

	
Orange	Black
	
Purple	Grey



It's tip of Nylon, high contact angle, and high friction coefficient EPDM, increases the user safety while decreasing the total weight.